



# Midline

## Venous Catheter Kit

**AMECATH**  
WE CARE


# Midline



## Venous Catheter Kits

**For non-pressure midline:** "Intended for vascular access peripherally in cephalic, basilica or the median cubital veins, basilica is the preferred cannulation site"

**For pressure resistance midline:** "Intended for vascular access cephalic, basilica or the median cubital veins, basilica is the peripherally in preferred cannulation site"  
"Identified by the purple color"



## "K" Configuration Midline Catheter Kit (Seldinger Technique) using IV cath cannula



### Kit Contents

K (Seldinger Technique)	
Component	Included
Midline Venous Catheter	✓
Over The Needle Catheter	✓
Straight End Nitinol Guide Wire	✓
Needless Connector (one for single lumen, and two for dual lumen)	✓

### Ordering Information

Midline "K" Kit Configuration

Single Lumen	Code	Fr Size	Length
	MLSLC-0310-K	3Fr	10 cm
	MLSLC-0312-K	3Fr	12 cm
	MLSLC-0320-K	3Fr	20 cm
	MLSLC-0325-K	3Fr	25 cm
	MLSLC-0410-K	4Fr	10 cm
	MLSLC-0412-K	4Fr	12 cm
	MLSLC-0420-K	4Fr	20 cm
	MLSLC-0425-K	4Fr	25 cm
	MLSLC-0512-K	5Fr	12 cm
	MLSLC-0520-K	5Fr	20 cm
	MLSLC-0525-K	5Fr	25 cm

### Ordering Information

Double Lumen	Code	Fr Size	Length
	MLDLC-0420-K	4Fr	20 cm
	MLDLC-0425-K	4Fr	25 cm
	MLDLC-0520-K	5Fr	20 cm
MLDLC-0525-K	5Fr	25 cm	

Pressure Resistant Midline "K" Kit Configuration

Pressure Single Lumen	Code	Fr Size	Length
	PRMLSLC-0310-K	3Fr	10 cm
	PRMLSLC-0312-K	3Fr	12 cm
	PRMLSLC-0320-K	3Fr	20 cm
	PRMLSLC-0325-K	3Fr	25 cm
	PRMLSLC-0410-K	4Fr	10 cm
	PRMLSLC-0412-K	4Fr	12 cm
	PRMLSLC-0420-K	4Fr	20 cm
	PRMLSLC-0425-K	4Fr	25 cm
	PRMLSLC-0512-K	5Fr	12 cm
PRMLSLC-0520-K	5Fr	20 cm	
PRMLSLC-0525-K	5Fr	25 cm	

Pressure Double	Code	Fr Size	Length
	PRMLDLC-0420-K	4Fr	20 cm
	PRMLDLC-0425-K	4Fr	25 cm
	PRMLDLC-0520-K	5Fr	20 cm
PRMLDLC-0525-K	5Fr	25 cm	

## "KON" Configuration Midline Catheter Kit for Seldinger Technique (using over needle peel away introducer)

### Kit Contents

KON (Direct Puncture)	
Component	Included
Midline Venous Catheter	✓
0.018" Nitinol Stiffener With Flexible End	✓
Over The Needle Peel Away Sheath Introducer	✓
Touhy Burst Adapter	✓
T-connector With Side Arm With Needless Connector	✓
10 CC Syringe	✓
Safety Scalpel	✓
Needless Connector (with dual lumen catheter only)	✓
Catheter Hub Securement Adhesive (Unifix Lock)	✓
Transparent Catheter Adhesive	✓



### Ordering Information

Midline "KON" Kit Configuration

Single Lumen	Code	Fr Size	Length
	MLSLC-0310-KON	3Fr	10 cm
	MLSLC-0312-KON	3Fr	12 cm
	MLSLC-0320-KON	3Fr	20 cm
	MLSLC-0325-KON	3Fr	25 cm
	MLSLC-0410-KON	4Fr	10 cm
	MLSLC-0412-KON	4Fr	12 cm
	MLSLC-0420-KON	4Fr	20 cm
	MLSLC-0425-KON	4Fr	25 cm
	MLSLC-0512-KON	5Fr	12 cm
	MLSLC-0520-KON	5Fr	20 cm
	MLSLC-0525-KON	5Fr	25 cm

Pressure Resistant Midline "KON" Kit Configuration

Single Lumen	Code	Fr Size	Length
	PRMLSLC-0310-KON	3Fr	10 cm
	PRMLSLC-0312-KON	3Fr	12 cm
	PRMLSLC-0320-KON	3Fr	20 cm
	PRMLSLC-0325-KON	3Fr	25 cm
	PRMLSLC-0410-KON	4Fr	10 cm
	PRMLSLC-0412-KON	4Fr	12 cm
	PRMLSLC-0420-KON	4Fr	20 cm
	PRMLSLC-0425-KON	4Fr	25 cm
	PRMLSLC-0512-KON	5Fr	12 cm
	PRMLSLC-0520-KON	5Fr	20 cm
	PRMLSLC-0525-KON	5Fr	25 cm

Double Lumen	Code	Fr Size	Length
	MLDLC-0420-KON	4Fr	20 cm
	MLDLC-0425-KON	4Fr	25 cm
	MLDLC-0520-KON	5Fr	20 cm
MLDLC-0525-KON	5Fr	25 cm	

Double Lumen	Code	Fr Size	Length
	PRMLDLC-0420-KON	4Fr	20 cm
	PRMLDLC-0425-KON	4Fr	25 cm
	PRMLDLC-0520-KON	5Fr	20 cm
PRMLDLC-0525-KON	5Fr	25 cm	

## "KOD" Configuration Midline Catheter Kit for modified Seldinger Technique (using Over Dilator Peel Away Introducer)

### Kit Contents

KOD (Modified Seldinger Technique)	
Component	Included
Midline Venous Catheter	√
0.018" Nitinol Stiffener With Flexible End	√
G21 X 7 Cm Introducer Needle With Echogenic Tip	√
Touhy Burst Adapter	√
T-connector With Side Arm With Needless Connector	√
10 Cc Syringe	√
Straight End 0.018"x 50 Cm Nitinol Guidewire	√
Over The Dilator Peel Away Sheath Introducer	√
Safety Scalpel	√
Needless Connector (with dual lumen catheter only)	√
Catheter Hub Securement Adhesive (Unifix Lock)	√
Transparent Catheter Adhesive	√

### Ordering Information

Midline "KOD" Kit Configuration

	Code	Fr Size	Length
Single Lumen	MLSLC-0310-KOD	3Fr	10 cm
	MLSLC-0312-KOD	3Fr	12 cm
	MLSLC-0320-KOD	3Fr	20 cm
	MLSLC-0325-KOD	3Fr	25 cm
	MLSLC-0410-KOD	4Fr	10 cm
	MLSLC-0412-KOD	4Fr	12 cm
	MLSLC-0420-KOD	4Fr	20 cm
	MLSLC-0425-KOD	4Fr	25 cm
	MLSLC-0512-KOD	5Fr	12 cm
	MLSLC-0520-KOD	5Fr	20 cm
	MLSLC-0525-KOD	5Fr	25 cm

	Code	Fr Size	Length
Double Lumen	MLDLC-0420-KOD	4Fr	20 cm
	MLDLC-0425-KOD	4Fr	25 cm
	MLDLC-0520-KOD	5Fr	20 cm
	MLDLC-0525-KOD	5Fr	25 cm


Pressure Resistant Midline "KOD" Kit Configuration

	Code	Fr Size	Length
Single Lumen	PRMLSLC-0310-KOD	3Fr	10 cm
	PRMLSLC-0312-KOD	3Fr	12 cm
	PRMLSLC-0320-KOD	3Fr	20 cm
	PRMLSLC-0325-KOD	3Fr	25 cm
	PRMLSLC-0410-KOD	4Fr	10 cm
	PRMLSLC-0412-KOD	4Fr	12 cm
	PRMLSLC-0420-KOD	4Fr	20 cm
	PRMLSLC-0425-KOD	4Fr	25 cm
	PRMLSLC-0512-KOD	5Fr	12 cm
	PRMLSLC-0520-KOD	5Fr	20 cm
	PRMLSLC-0525-KOD	5Fr	25 cm

	Code	Fr Size	Length
Double Lumen	PRMLDLC-0420-KOD	4Fr	20 cm
	PRMLDLC-0425-KOD	4Fr	25 cm
	PRMLDLC-0520-KOD	5Fr	20 cm
	PRMLDLC-0525-KOD	5Fr	25 cm

## KOD Continued


G21 X 7 CM INTRODUCER  
NEEDLE WITH ECHOGENIC TIP



Straight end guide wire with  
Roller dispenser

PEEL-AWAY OVER DILATOR  
SHEATH INTRODUCER





NEEDLESS CONNECTOR

T-CONNECTOR

TOUHY CONNECTOR

NITINOL STIFFENER



2797

Manufactured By:  
**AMECO MEDICAL INDUSTRIES**

Industrial Zone B4 Plot 119 east,  
10th of Ramadan city - Egypt  
Tel: +20.554.501321/ 2  
Fax:+20.554.501224  
Email: support@amecathgroup.com  
Website: www.amecathgroup.com

**AMECATH**  
WE CARE

[www.amecathgroup.com](http://www.amecathgroup.com)



Authorized European Representative  
**Obelis s.a**

Bd. Général Wahis 53  
1030 Brussels, BELGIUM  
Tel: +32.2. 732.59.54  
Fax: +32.2.732.60.03  
Email: mail@obelis.net  
Website: www.obelis.net