



AMECATH
WE CARE

www.amecathgroup.com



2797

Manufactured By:
AMECO MEDICAL INDUSTRIES

Industrial Zone B4 Plot 119 east,
10th of Ramadan city - Egypt
Tel: +20.554.501321/ 2
Fax: +20.554.501224
Email: support@amecathgroup.com
Website: www.amecathgroup.com

Authorized European Representative
Obelis s.a

Bd. Général Wahis 53
1030 Brussels, BELGIUM
Tel: +32.2. 732.59.54
Fax: +32.2.732.60.03
Email: mail@obelis.net
Website: www.obelis.net

PICC

Peripheral Inserted
Central Venous Catheter Kit

AMECATH
WE CARE

PICC

Peripheral Inserted
Central Venous Catheter Kit



For Non-Pressure PICC: "Intended for peripheral vascular access"

For Pressure PICC: "Intended for peripheral vascular access & Power Injection of contrast media"



"KON" Configuration

Direct Puncture Technique Using Over The Needle Peel-Away Sheath Introducer



Kit Contents

Power resistant PICC

KON (DIRECT PUNCTURE)	
COMPONENT	INCLUDED
Peripherally Inserted Central Venous Catheter	√
Touhy-burst Adaptor	√
T-connector With Sidearm With Needleless Connector	√
0.018" Nitinol Stiffener With Flexible Tip	√
Over The Needle Peel Away Sheath Introducer	√
Safety Scalpel	√
10 Cc Syringe Lock	√
Needleless Connector (Valve) – one for dual lumen and two for triple lumen	√
Catheter Hub Securement Adhesive (Unifix Lock)	√
Transparent Catheter Adhesives	√
Measuring Tape	√

Ordering Information

	Code	Fr Size	Length
Single Lumen	SLPICC-0340-KON	3	40
	SLPICC-0350-KON	3	50
	SLPICC-0355-KON	3	55
	SLPICC-0360-KON	3	60
	SLPICC-0440-KON	4	40
	SLPICC-0450-KON	4	50
	SLPICC-0455-KON	4	55
	SLPICC-0460-KON	4	60
	SLPICC-0540-KON	5	40
	SLPICC-0550-KON	5	50
SLPICC-0555-KON	5	55	
SLPICC-0560-KON	5	60	

	Code	Fr Size	Length
Double Lumen	DLPICC-0440-KON	4	40
	DLPICC-0450-KON	4	50
	DLPICC-0455-KON	4	55
	DLPICC-0460-KON	4	60
	DLPICC-0540-KON	5	40
	DLPICC-0550-KON	5	50
	DLPICC-0555-KON	5	55
	DLPICC-0560-KON	5	60
	DLPICC-0640-KON	6	40
	DLPICC-0650-KON	6	50
DLPICC-0655-KON	6	55	
DLPICC-0660-KON	6	60	

	Code	Fr Size	Length
Triple Lumen	TLPICC-0640-KON	6	40
	TLPICC-0650-KON	6	50
	TLPICC-0655-KON	6	55
	TLPICC-0660-KON	6	60


	Code	Fr Size	Length
Pressure Single Lumen	PRSLPICC-0355-KON	3	55
	PRSLPICC-0455-KON	4	55
	PRSLPICC-0555-KON	5	55

	Code	Fr Size	Length
Pressure Double Lumen	PRDLPICC-0555-KON	5	55
	PRDLPICC-0655-KON	6	55

	Code	Fr Size	Length
Pressure Triple Lumen	PRTLPICC-0655-KON	6	55

"KOD" Configuration

Modified Seldinger Technique Using Over The Dilator Peel-Away Sheath Introducer



KOD (MODIFIED SELDINGER TECHNIQUE)	
COMPONENT	INCLUDED
Peripherally Inserted Central Venous Catheter	√
Touhy-burst Adaptor	√
T-connector With Sidearm With Needleless Connector	√
0.018" Nitinol Stiffener With Flexible Tip	√
G21x 7 Cm Introducer Needle With Echogenic Tip	√
10 Cc Syringe Lock	√
Straight End 0.018" X 50 Cm Nitinol Core Guidewires	√
OVER THE DILATOR PEEL AWAY SHEATH INTRODUCER	√
Safety Scalpel	√
NEEDLELESS CONNECTOR (Valve) – One For Dual Lumen And Two For Triple Lumen	√
Catheter Hub Securement Adhesive (Unifix Lock)	√
Transparent Catheter Adhesive	√
Measuring Tape	√

Ordering Information

	Code	Fr Size	Length
Single Lumen	SLPICC-0340-KOD	3	40
	SLPICC-0340-KOD+	3	40
	SLPICC-0350-KOD	3	50
	SLPICC-0350-KOD+	3	50
	SLPICC-0355-KOD	3	55
	SLPICC-0355-KOD+	3	55
	SLPICC-0360-KOD	3	60
	SLPICC-0360-KOD+	3	60
	SLPICC-0440-KOD	4	40
	SLPICC-0440-KOD+	4	40
	SLPICC-0450-KOD	4	50
	SLPICC-0450-KOD+	4	50
	SLPICC-0455-KOD	4	55
	SLPICC-0455-KOD+	4	55
	SLPICC-0460-KOD	4	60
	SLPICC-0460-KOD+	4	60
SLPICC-0540-KOD	5	40	
SLPICC-0550-KOD	5	50	
SLPICC-0555-KOD	5	55	
SLPICC-0560-KOD	5	60	

	Code	Fr Size	Length
Double Lumen	DLPICC-0440-KOD	4	40
	DLPICC-0440-KOD+	4	40
	DLPICC-0450-KOD	4	50
	DLPICC-0450-KOD+	4	50
	DLPICC-0455-KOD	4	55
	DLPICC-0455-KOD+	4	55
	DLPICC-0460-KOD	4	60
	DLPICC-0460-KOD+	4	60
	DLPICC-0540-KOD	5	40
	DLPICC-0550-KOD	5	50
	DLPICC-0555-KOD	5	55
	DLPICC-0560-KOD	5	60
	DLPICC-0640-KOD	6	40
	DLPICC-0650-KOD	6	50
DLPICC-0655-KOD	6	55	
DLPICC-0660-KOD	6	60	

	Code	Fr Size	Length
Triple Lumen	TLPICC-0640-KOD	6	40
	TLPICC-0650-KOD	6	50
	TLPICC-0655-KOD	6	55
	TLPICC-0660-KOD	6	60

Add "+" Sign at the end of some codes represents configurations in which the peel away over dilator sheath is over-sized by 0.5Fr.

	Code	Fr Size	Length
Pressure Single Lumen	PRSLPICC-0355-KOD	3	55
	PRSLPICC-0355-KOD+	3	55
	PRSLPICC-0455-KOD	4	55
	PRSLPICC-0455-KOD+	4	55
	PRSLPICC-0555-KOD	5	55

	Code	Fr Size	Length
Pressure Double Lumen	PRDLPICC-0555-KOD	5	55
	PRDLPICC-0655-KOD	6	55

	Code	Fr Size	Length
Pressure Triple Lumen	PRTL PICC-0655-KOD	6	55