



AMECATH
WE CARE

www.amecathgroup.com

Pediatric Central Venous CATHETER KIT



2797

Manufactured By:
AMECO MEDICAL INDUSTRIES

Industrial Zone B4 Plot 119 east,
10th of Ramadan city - Egypt
Tel: +20.554.501321/ 2
Fax: +20.554.501224
Email: support@amecathgroup.com
Website: www.amecathgroup.com

Authorized European Representative
Obelis s.a

Bd. Général Wahis 53
1030 Brussels, BELGIUM
Tel: +32.2. 732.59.54
Fax: +32.2.732.60.03
Email: mail@obelis.net
Website: www.obelis.net

AMECATH
WE CARE


NEOLINE



" KP Configuration " Over The Wire Technique

" Mainly intended for vascular access by the jugular or subclavian route "

Mini Female Luer
Extra Comfort





Extension Line
Integral Extension Lines
Provides less
manipulation at the
insertion Site

**Mini Hub
With Fixation Wings**
Designed specifically for
use with Neonates and
Pediatric Patients



PUR catheter
marked every 1 cm in
order to facilitate an
easy, precise and safe
positioning

Catheter Tip
Tapered For Smooth
Insertion and Soft to
Protect the Vessel Wall



Kit Contents

| KP (over the wire Technique) | |
|--------------------------------------|----------|
| Component | Included |
| Polyurethane Central Venous Catheter | ✓ |
| Guide Wire | ✓ |
| Vessel Dilator | ✓ |
| Syringe | ✓ |
| Introducer Needle | ✓ |
| Scalpel | ✓ |
| Movable Fixation Wing | ✓ |
| Injection caps | ✓ |
| Catheter Over The Needle | Optional |
| Neutral Valves | Optional |
| Guiding Y Connector | Optional |



**Guiding Y
Connector**



**Catheter over
the needle**



Neutral Valve

- Add " O " To the end of The Code in case of Catheter over the Needle instead of the Standard Introducer Needle .
- Add " V " To end of The Code in case of Neutral Valves .
- Add " GY " To end of The Code in case of guiding Y Connector .
- Add " A " to the End of the Code in case of Guide Wire with Advancer.

Ordering information

| Single Lumen | Code | Fr Size | Length |
|---------------|---------------|---------|--------|
| | CSLC- 0205-KP | 2 | 5 |
| | CSLC- 0208-KP | 2 | 8 |
| | CSLC- 0305-KP | 3 | 5 |
| | CSLC- 0308-KP | 3 | 8 |
| | CSLC- 0310-KP | 3 | 10 |
| | CSLC- 0313-KP | 3 | 13 |
| | CSLC- 0315-KP | 3 | 15 |
| | CSLC- 0317-KP | 3 | 17 |
| | CSLC- 0405-KP | 4 | 5 |
| | CSLC- 0406-KP | 4 | 6 |
| | CSLC- 0408-KP | 4 | 8 |
| | CSLC- 0410-KP | 4 | 10 |
| | CSLC- 0412-KP | 4 | 12 |
| | CSLC- 0413-KP | 4 | 13 |
| | CSLC- 0415-KP | 4 | 15 |
| | CSLC- 0417-KP | 4 | 17 |
| CSLC- 0420-KP | 4 | 20 | |
| CSLC- 0506-KP | 5 | 6 | |
| CSLC- 0508-KP | 5 | 8 | |
| CSLC- 0510-KP | 5 | 10 | |
| CSLC- 0513-KP | 5 | 13 | |
| CSLC- 0515-KP | 5 | 15 | |
| CSLC- 0517-KP | 5 | 17 | |
| CSLC- 0520-KP | 5 | 20 | |

| Double Lumens | Code | Fr Size | Length |
|-----------------|-----------------|---------|--------|
| | CDLC- 0304-KP | 3 | 4 |
| | CDLC- 0305-KP | 3 | 5 |
| | CDLC- 0308-KP | 3 | 8 |
| | CDLC- 0310-KP | 3 | 10 |
| | CDLC- 0313-KP | 3 | 13 |
| | CDLC- 0405-KP | 4 | 5 |
| | CDLC- 0406-KP | 4 | 6 |
| | CDLC- 0408-KP | 4 | 8 |
| | CDLC- 0410-KP | 4 | 10 |
| | CDLC- 0413-KP | 4 | 13 |
| | CDLC- 0415-KP | 4 | 15 |
| | CDLC- 04.508-KP | 4.5 | 8 |
| | CDLC- 04.510-KP | 4.5 | 10 |
| | CDLC- 0505-KP | 5 | 5 |
| | CDLC- 0506-KP | 5 | 6 |
| | CDLC- 0508-KP | 5 | 8 |
| CDLC- 0510-KP | 5 | 10 | |
| CDLC- 0513-KP | 5 | 13 | |
| CDLC- 0515-KP | 5 | 15 | |
| CDLC- 05.510-KP | 5.5 | 10 | |
| CDLC- 05.515-KP | 5.5 | 15 | |

| Triple Lumens | Code | Fr Size | Length |
|-----------------|---------------|---------|--------|
| | CTLC- 0405-KP | 4 | 5 |
| | CTLC- 0406-KP | 4 | 6 |
| | CTLC- 0408-KP | 4 | 8 |
| | CTLC- 0410-KP | 4 | 10 |
| | CTLC- 0413-KP | 4 | 13 |
| | CTLC- 0430-KP | 4 | 30 |
| | CTLC- 0505-KP | 5 | 5 |
| | CTLC- 0508-KP | 5 | 8 |
| | CTLC- 0510-KP | 5 | 10 |
| | CTLC- 0513-KP | 5 | 13 |
| | CTLC- 0515-KP | 5 | 15 |
| | CTLC- 0520-KP | 5 | 20 |
| CTLC- 0530-KP | 5 | 30 | |
| CTLC- 05.510-KP | 5.5 | 10 | |
| CTLC- 05.515-KP | 5.5 | 10 | |

| Quadra Lumen | Code | Fr Size | Length |
|--------------|--------------|---------|--------|
| | CQLC-0608-KP | 6 | 8 |
| | CQLC-0613-KP | 6 | 13 |
| | CQLC-0615-KP | 6 | 15 |
| CQLC-0620-KP | 6 | 20 | |

" KP* " Configuration 'Peripherally inserted catheter using the Peel-able cannula (Direct Cannulation Technique)

The peel able cannula enables the catheter to be safely advanced through plastic. Additional security is provided by a small primary wing to help minimize the risk of catheter movement and kinking.

Kit Contents

| KP* (Peripheral CVC) | |
|--|----------|
| Component | Included |
| Polyurethane Central Venous Catheter | ✓ |
| Peel-away over the Needle | ✓ |
| Injection caps | ✓ |
| Measuring tape | ✓ |
| T-connector with Stiffener Wire & Caps | ✓ |

Ordering information

| Single Lumen | Code | Fr Size | Length |
|----------------|----------------|---------|--------|
| | CSLC- 0205-KP* | 2 | 5 |
| | CSLC- 0208-KP* | 2 | 8 |
| | CSLC- 0215-KP* | 2 | 15 |
| | CSLC- 0230-KP* | 2 | 30 |
| | CSLC- 0305-KP* | 3 | 5 |
| | CSLC- 0308-KP* | 3 | 8 |
| | CSLC- 0313-KP* | 3 | 13 |
| CSLC- 0315-KP* | 3 | 15 | |
| CSLC- 0330-KP* | 3 | 30 | |



" KP " Configuration 'Peripherally inserted catheter using the Peel-away over Dilator (modified Seldinger Technique)**
The Peel-able Over Dilator sheath enables the catheter to be safely advanced.



Kit Contents

| KP** (Peripheral CVC) | |
|--|----------|
| Component | Included |
| Polyurethane Central Venous Catheter | ✓ |
| Guide Wire with Dispenser | ✓ |
| Peel-away over Dilator | ✓ |
| Syringe | ✓ |
| Introducer Needle | ✓ |
| Measuring Tape | ✓ |
| Injection caps | ✓ |
| T-connector with Stiffener Wire & Caps | ✓ |

Ordering information

| | Code | Fr Size | Length |
|-------------------------------|-----------------|---------|--------|
| Double Lumens Single Lumen | CSLC- 0305-KP** | 3 | 5 |
| | CSLC- 0308-KP** | 3 | 8 |
| | CSLC- 0313-KP** | 3 | 13 |
| | CSLC- 0315-KP** | 3 | 15 |
| | CSLC- 0330-KP** | 3 | 30 |
| | CDLC- 0305-KP** | 3 | 5 |
| | CDLC- 0308-KP** | 3 | 8 |
| | CDLC- 0313-KP** | 3 | 13 |
| | CDLC- 0315-KP** | 3 | 15 |
| | CDLC- 0330-KP** | 3 | 30 |

" KU" Configuration Umbilical Catheterization - Cut Down Technique

The safest choice for short-term vascular access in neonates

Kit Contents

| Umbilical | |
|--|----------|
| Component | Included |
| Polyurethane Catheter | ✓ |
| 3-Way Stopcock | ✓ |
| T-connector with Stiffener Wire & Caps (2F & 3F) ONLY | ✓ |
| Injection caps | ✓ |

Ordering information

| | Code | Fr Size | Length |
|---------------|---------------|---------|--------|
| Single Lumen | CSLC- 0240-KU | 2 | 40 |
| | CSLC- 0340-KU | 3 | 40 |
| | CSLC-0440-KU | 4 | 40 |
| | CSLC- 0540-KU | 5 | 40 |
| | CSLC- 0840-KU | 8 | 40 |
| Double Lumens | CDLC-0325-KU | 3 | 25 |
| | CDLC-0340-KU | 3 | 40 |
| | CDLC-0425-KU | 4 | 25 |
| | CDLC-0440-KU | 4 | 40 |
| | CDLC-0525-KU | 5 | 25 |
| Triple Lumens | CDLC-0540-KU | 5 | 40 |
| | CTLC-0425-KU | 4 | 25 |

