



Short Term

Hemodialysis Catheters Kit

AMECATH
WE CARE

1. Desilet Kit (DK)



" Device placed in the vascular system by the Seldinger technique for introducing cardiovascular or silicone catheters "

Kit Contents


- Polyurethane Single lumen Haemodialysis catheter.
- J End Nitinol guide wire.
- Vessel dilator
- Syringe
- Introducer needle, echogenic
- Scalpel
- Injection cap
- Towel



Ordering Information

Fr. Size	Length	Desilet Kit Code
10	13 Cm	SLC-1013-DK
10	15 Cm	SLC-1015-DK
10	18 Cm	SLC-1018-DK
10	20 Cm	SLC-1020-DK

2. Tapered Tip Double Lumens Catheters Kit (DLC)



Guiding Syringe



Soft Tapered Tip

- Atraumatic Tip , reduces Potential vessel trauma



Guiding Y Connector

Kit Contents

Standard Kits • Tapered Tip Dual	Standard Kits with Guiding syringe
Lumen Haemodialysis • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilators • Syringe • Introducer Needle, Echogenic • Scalpel • Injection Cap • Movable Fixation Wings (with pre-curved catheter) • Guiding Y Connector	• Tapered Tip Dual Lumen Haemodialysis • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilators • Guiding Syringe • Introducer Needle, Echogenic • Scalpel • Injection Cap • Movable Fixation Wings (with pre-curved catheter)

Subclavian Catheter




Subclavian Catheter



Jugular Pre Curved Catheter



Kit With Guiding Syringe



Ordering Information

Fr. Size	Length	Subclavian Kit Code	Jugular Kit Code	Pre-curved Kit Code
06	08Cm	DLC-0608-KP	XX	XX
07	10Cm	DLC-0710-KP	DLC-0710-KJP	XX
07	13Cm	DLC-0713-KP	DLC-0713-KJP	XX
08	10Cm	DLC-0810-KP	DLC-0810-KJP	XX
08	12Cm	DLC-0812-KP	DLC-0812-KJP	XX
08	13Cm	DLC-0813-KP	DLC-0813-KJP	XX
09	10Cm	DLC-0910-KP	DLC-0910-KJP	XX
09	12Cm	DLC-0912-KP	DLC-0912-KJP	XX
09	13Cm	DLC-0913-KP	DLC-0913-KJP	XX
09	13Cm	DLC-0913-K	DLC-0913-KJ	DLC-0913-KPC
09	15Cm	DLC-0915-K	DLC-0915-KJ	DLC-0915-KPC
10	10Cm	DLC-1010-K	DLC-1010-KJ	XX
10	13Cm	DLC-1013-K	DLC-1013-KJ	DLC-1013-KPC
10	15Cm	DLC-1015-K	DLC-1015-KJ	DLC-1015-KPC
12	12Cm	DLC-1212-K	DLC-1212-KJ	DLC-1212-KPC
12	13Cm	DLC-1213-K	DLC-1213-KJ	DLC-1213-KPC
12	15Cm	DLC-1215-K	DLC-1215-KJ	DLC-1215-KPC
12	16Cm	DLC-1216-K	DLC-1216-KJ	DLC-1216-KPC
12	17Cm	DLC-1217-K	DLC-1217-KJ	DLC-1217-KPC
12	20Cm	DLC-1220-K	DLC-1220-KJ	DLC-1220-KPC
12	24Cm	DLC-1224-K	DLC-1224-KJ	DLC-1224-KPC
14	15Cm	DLC-1415-K	DLC-1415-KJ	DLC-1415-KPC
14	17Cm	DLC-1417-K	DLC-1417-KJ	DLC-1417-KPC
14	20Cm	DLC-1420-K	DLC-1420-KJ	DLC-1420-KPC
16	15Cm	DLC-1615-K	DLC-1615-KJ	XX
16	17Cm	DLC-1617-K	DLC-1617-KJ	XX
16	20Cm	DLC-1620-K	DLC-1620-KJ	XX

Add " GS " To the End of the Code for the Kit with a Guiding Syringe .

3. Step Tip Double Lumens Catheters Kit (SDLC)



Step Tip Double Lumens Catheter



Step Tip Design

- For maximum flow
- Decrease Recirculation

Kit Contents

Standard Kits	Standard Kits with Guiding syringe
<ul style="list-style-type: none"> • Step Tip Soft Dual Lumen Haemodialysis • Stiffener with Haemostatic Valve • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilator • Vessel Dilator with hydrophilic coating • Syringe • Guiding Y Connector • Introducer Needle, Echogenic • Scalpel • Injection Cap • Unifix Catheter Tube Fixation Adhesive • Transparent Catheter Fixation Adhesive • Movable Fixation Wings 	<ul style="list-style-type: none"> • Step Tip Soft Dual Lumen Haemodialysis • Stiffener with Haemostatic Valve • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilator • Vessel Dilator with hydrophilic coating • Guiding Syringe • Introducer Needle, Echogenic • Scalpel • Injection Cap • Unifix Catheter Tube Fixation Adhesive • Transparent Catheter Fixation Adhesive • Movable Fixation Wings

Ordering Information

Fr. Size	Length	Subclavian Kit Code	Jugular Kit Code	Pre-curved Kit Code
12	12Cm	SDLC-1212-K	SDLC-1212-KJ	XX
12	15Cm	SDLC-1215-K	SDLC-1215-KJ	SDLC-1215-KPC
12	17Cm	SDLC-1217-K	SDLC-1217-KJ	SDLC-1217-KPC
12	20Cm	SDLC-1220-K	SDLC-1220-KJ	SDLC-1220-KPC
14	25Cm	SDLC-1225-K	SDLC-1225-KJ	SDLC-1225-KPC
14	15Cm	SDLC-1415-K	SDLC-1415-KJ	SDLC-1415-KPC
14	17Cm	SDLC-1417-K	SDLC-1417-KJ	SDLC-1417-KPC
14	20Cm	SDLC-1420-K	SDLC-1420-KJ	SDLC-1420-KPC
14	25Cm	SDLC-1425-K	SDLC-1425-KJ	SDLC-1425-KPC
14	30Cm	SDLC-1430-K	SDLC-1430-KJ	SDLC-1430-KPC
14	35Cm	SDLC-1435-K	SDLC-1435-KJ	SDLC-1435-KPC

Add " GS " To the End of the Code for the Kit with a Guiding Syringe.



Kit With Guiding Syringe

4. SMART DUAL LUMENS SHORT TERM HAEMODIALYSIS CATHETER Kit (DLC-SM)



Smart Grooves Design Prevents the side Holes from being Stick to the vessel wall During blood Suction.

Kit Contents

Standard Kits	Standard Kits with Guiding Syringe
<ul style="list-style-type: none"> • Dual lumen Smart Haemodialysis Catheter • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilators • Syringe • Guiding Y Connector • Introducer Needle, Echogenic • Scalpel • Injection Cap • Movable fixation Wing (with pre-curved catheter) 	<ul style="list-style-type: none"> • Dual lumen Smart Haemodialysis Catheter • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilators • Guiding Syringe • Introducer Needle, Echogenic • Scalpel • Injection Cap • Movable fixation Wing (with pre-curved catheter)

Ordering Information

Fr. Size	Length	Subclavian Kit Code	Jugular Kit Code	Pre-curved Kit Code
12	12Cm	DLC-1212-KSM	DLC-1212-KJSM	DLC-1212-KPCSM
12	15Cm	DLC-1215-KSM	DLC-1215-KJSM	DLC-1215-KPCSM
12	17Cm	DLC-1217-KSM	DLC-1217-KJSM	DLC-1217-KPCSM
12	20Cm	DLC-1220-KSM	DLC-1220-KJSM	DLC-1220-KPCSM
14	15Cm	DLC-1415-KSM	DLC-1415-KJSM	DLC-1415-KPCSM
14	17Cm	DLC-1417-KSM	DLC-1417-KJSM	DLC-1417-KPCSM
14	20Cm	DLC-1420-KSM	DLC-1420-KJSM	DLC-1420-KPCSM

5. UniQaih Dual Lumens Catheters Kit (UN-DLC)

- **Soft material:** comfortable for the patient & easily inserted.
- **Smart groove:** protect side holes from sticking to vessel wall.
- **Tubular design:** No turbulence, higher flow
- **No connection:** Prevent turbulence, more comfort
- **No clamp:** saved integrity of extensions.
- **Hemostatic valves:** prevent blood leakage all over the procedure.



Kit Contents

Standard Kits	Standard Kits with Guiding Syringe
<ul style="list-style-type: none"> • Dual Lumen Haemodialysis Catheter • Stiffener with Haemostatic Valve • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilator • Vessel Dilator with hydrophilic coating • Syringe • Guiding Y Connector • Introducer Needle, Echogenic • Scalpel • Injection Cap • Unifix Catheter Tube Fixation Adhesive • Transparent Catheter Fixation Adhesive • Movable fixation wings 	<ul style="list-style-type: none"> • Dual Lumen Haemodialysis Catheter • Stiffener with Haemostatic Valve • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilator • Vessel Dilator with hydrophilic coating • Guiding Syringe • Introducer Needle, Echogenic • Scalpel • Injection Cap • Unifix Catheter Tube Fixation Adhesive • Transparent Catheter Fixation Adhesive • Movable Fixation wings.

Ordering Information

Fr Size	Length	Subclavian Kit Code
12	15Cm	UN-DLC-1215-K
12	17Cm	UN-DLC-1217-K
12	20Cm	UN-DLC-1220-K
12	25Cm	UN-DLC-1225-K
14	15Cm	UN-DLC-1415-K
14	17Cm	UN-DLC-1417-K
14	20Cm	UN-DLC-1420-K
14	25Cm	UN-DLC-1425-K



Smart Grooves Design Prevents the side Holes from being Stick to the vessel wall During blood Suction.

Add " GS " To the End of the Code for the kit with a Guiding Syringe.

6. Triple Lumens Catheters Kit (TLC)

Third lumen gives greater flexibility for fluid administration, vasoactive drug therapy, nutrition and blood sampling



Kit Contents

Standard Kits	Standard Kits with Guiding Syringe
<ul style="list-style-type: none"> • Triple Lumen Haemodialysis Catheter • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilators • Syringe • Guiding Y Connector • Introducer Needle, Echogenic • Scalpel • Injection Cap • Movable fixation Wing (with pre-curved catheter) 	<ul style="list-style-type: none"> • Triple Lumen Haemodialysis Catheter • J End Nitinol Guide Wire with Scaled Dispenser • Vessel Dilators • Guiding Syringe • Introducer Needle, Echogenic • Scalpel • Injection Cap • Movable fixation Wing (with pre-curved catheter)


Ordering Information

Fr. Size	Length	Subclavian Kit Code	Jugular Kit Code	Pre-curved Kit Code
13	15Cm	TLC-1215-K	TLC-1215-KJ	TLC-1215-KPC
13	15Cm	TLC-1315-K	TLC-1315-KJ	TLC-1315-KPC
13	17Cm	TLC-1317-K	TLC-1317-KJ	TLC-1317- KPC
13	20Cm	TLC-1320-K	TLC-1320-KJ	TLC-1320- KPC
13	24Cm	TLC-1324-K	TLC-1324-KJ	TLC-1324-KPC
13	25Cm	TLC-1325-K	TLC-1325-KJ	TLC-1325-KPC
14	15Cm	TLC-1415-K	TLC-1415-KJ	TLC-1415- KPC
14	17Cm	TLC-1417-K	TLC-1417-KJ	TLC-1417- KPC
14	20Cm	TLC-1420-K	TLC-1420-KJ	TLC-1420- KPC

Add " GS " To the End of the Code for the Kit with a Guiding Syringe.

7. SMART Triple Lumens Catheter Kit (TLC-SM)

Third lumen gives greater flexibility for fluid administration, vasoactive drug therapy, nutrition and blood sampling



Smart Grooves Design Prevents the side Holes from being Stick to the vessel wall During blood Suction.

Kit Contents

Standard Kits	Standard Kits with Guiding Syringe
<ul style="list-style-type: none"> • Triple Lumen Haemodialysis Catheter • J End Nitinol Guide Wire with Scaled <p>Dispenser</p> <ul style="list-style-type: none"> • Vessel Dilators • Syringe • Guiding Y Connector • Introducer Needle, Echogenic • Scalpel • Injection Cap • Movable fixation Wing (with Pre-curved catheter) 	<ul style="list-style-type: none"> • Triple Lumen Haemodialysis Catheter • J End Nitinol Guide Wire with Scaled <p>Dispenser</p> <ul style="list-style-type: none"> • Vessel Dilators • Guiding Syringe • Introducer Needle, Echogenic • Scalpel • Injection Cap • Movable fixation Wing (with Pre-curved catheter)

Ordering Information

Fr. Size	Length	Subclavian Kit Code	Jugular Kit Code	Pre-curved Kit Code
13	15Cm	TLC-1315-KSM	TLC-1315-KJSM	TLC-1315-KPCSM
13	17Cm	TLC-1317- KSM	TLC-1317- KJSM	TLC-1317- KPCSM
13	20Cm	TLC-1320- KSM	TLC-1320- KJSM	TLC-1320- KPCSM
14	15Cm	TLC-1415- KSM	TLC-1415- KJSM	TLC-1415- KPCSM
14	17Cm	TLC-1417- KSM	TLC-1417- KJSM	TLC-1417- KPCSM
14	20Cm	TLC-1420- KSM	TLC-1420- KJSM	TLC-1420- KPCSM

Add " GS " To the End of the Code for the Kit with a Guiding Syringe.



2797

Manufactured By:

AMECO MEDICAL INDUSTRIES

Industrial Zone B4 Plot 119 east,
10th of Ramadan city - Egypt

Tel:+20.554.501321/2

Fax:+20.554.501224

Email:support@amecathgroup.com

Website: www.amecathgroup.com

AMECATH
WE CARE

www.amecathgroup.com

Authorized European Representative

Obelis S.A.

Bd. Général Wahis 53

1030 Brussels, BELGIUM

Tel: +32.2. 732.59.54

Fax: +32.2.732.60.03

Email: mail@obelis.net

Website: www.obelis.net